

INC-8 biological incubator is indicated to check the correct operation of ethylene oxide sterilizers and autoclaves.

It works by means of a simple procedure: a biological indicator is introduced in a normal sterilization cycle in the autoclave or sterilizer, and when finished, it is introduced into the INC-8 device.

This operation must be performed regularly to ensure successful sterilization.



FEATURES

Incubator for biological control.

Simple handling.

For ethylene oxide sterilizers and autoclaves.

Compact structure.

Temperature from 37°C to 57°C.

Dimensions: 180x175x80mm.

Nominal power: 50W.

Incubation of 8 vials simultaneously.